NCTF® 135
Acide hyaluronique
+ solution polyrévitalisante

NCTF® 135HA
Acide hyaluronique hautement concentré
+ solution polyrévitalisante

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DESCRIPTION

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CONDITIONS DE STOCKAGE / CONSERVATION

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INSTRUCTIONS FOR USE

DESCRIPTION

5 x 3ml vials
5 syringes of 3ml
5 x 18G withdrawal needles
5 x 30G treatment needles
5 x 32G treatment needles

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5 Nucleic acids

Deoxyadenosine
Deoxycytidine
Deoxyguanosine
Deoxythymidine

5-Methyl-2'-deoxycytidine

24 Amino Acids

Acid α -aminobutyric
Alanine
Arginine
Asparagine
Aspartic acid
Cystine
Glutamine
Glutamic acid
Glycine
Histidine
Hydroxyproline
Isoleucine
Leucine
Lysine
Methionine
Ornithine
Phenylalanine
Proline
Serine
Taurine
Threonine
Tryptophan
Tyrosine
Valine

6 Coenzymes

TPP (Carboxylase)
CoA (Coenzyme A)
FAD (Flavin adenine dinucleotide)
NAD (Nicotinamide adenine dinucleotide)
NADP (Nicotinamide adenine dinucleotide phosphate)
UTP (Uridine triphosphate)
Other
Polysorbitate 80
Acid glucuronic
Acid glucuronic lactone
Glucosamine
Dextrose anhydrous

INDICATIONS

NCTF® 135: Anti-ageing mesotherapy product indicated for revitalization, hydration of tired or dull skin and treatment of superficial wrinkles.

NCTF® 135HA: Anti-ageing mesotherapy product indicated for intense revitalization and hydration of tired or dull skin, treatment of wrinkles and redensification of mature or slack skin.

NCTF® 135 and NCTF® 135HA are indicated to be injected in the epidermis or dermis in the following areas:

- Face (including the peri-orbital area)
- Neck
- Décolleté area
- Back of the hands
- Internal face of the arms

CONTRAINDICATIONS

No inject in blood vessels. Do not use in:

- patients with a skin alteration or disease of any kind in the treated area.
- patients with a history of auto-immune disease or under immunotherapy.
- patients with hypersensitivity or a known allergy to any of the constituents.
- pregnant or breast-feeding women.
- children.

PRECAUTIONS FOR USE

There are no clinical data available in terms of tolerance with regard to the injection of NCTF® 135 and NCTF® 135HA in an area already treated with another "permanent" filling product. Precaution should be taken in such cases, although NCTF® 135 and NCTF® 135HA are not injected at the same level.

Patients receiving an anticoagulant treatment or a platelet aggregation inhibitor (for example, aspirin) must be informed of the increased risk of bruising and mild bleeding at the injection site.

The patient should be given the recommendation not to use sunscreens or tanning beds following the injection and to avoid prolonged exposure to the sun, UV rays, freezing temperatures as well as use of saunas or steam rooms for at least 48 hours after the injection.

Training in mesotherapy techniques is recommended for the use of NCTF® 135 and NCTF® 135HA.

A local anaesthetic (patch or cream) can be applied prior to the injection. Please note that anaesthetic products can cause local redness or hypersensitivity.

Do not inject NCTF® 135 or NCTF® 135HA into beauty spots (navel) and in lesions cutanées de type acné actif.

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There is a known incompatibility between hyaluronic acid and quaternary ammonium salts such as benzalkonium chloride.

NCTF® 135 and NCTF® 135HA must not be injected into contact with such products or with medical/surgical equipment treated with this type of product.

UNDESIRABLE EFFECTS

The practitioner must inform the patient that there are possible undesirable effects related to the implantation of this device, such as:

• Redness or mild local inflammation that usually disappears within 24 or 48 hours.

• Pain or discomfort at the injection site that stops rapidly and spontaneously after the injection.

• Light oedema and small haematomas may occur, but disappear within 48 hours.

• Pain or discomfort at the injection site that stops rapidly and spontaneously after the injection.

• Transient pain at the injection site.

• Slight bleeding or bruising following the injection and to avoid prolonged exposure to the sun, UV rays, freezing temperatures as well as use of saunas or steam rooms for at least 48 hours after the injection.

• Training in mesotherapy techniques is recommended for the use of NCTF® 135 and NCTF® 135HA.

A local anaesthetic (patch or cream) can be applied prior to the injection. Please note that anaesthetic products can cause local redness or hypersensitivity.

Do not inject NCTF® 135 or NCTF® 135HA into beauty spots (navel) and in lesions cutanées de type acné actif.

INSTRUCTIONS FOR USE

Before starting treatment, the patient must be informed of the indications of the device, its contraindications, its incompatibilities and its possible undesirable effects.

All asepsis and hygiene precautions must be taken to prepare the site, the material and the patient prior to the injection.

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WARNINGS

- Check the expiry date on the label.
- Do not reuse.
- Do not reuse the product.
- At use, the needles must be thrown away in a suitable safe container.
- After use, the needles must be disposed in a suitable safe container.
- Do not use if the packaging is open or damaged

STORAGE / CONSERVATION CONDITIONS

Store NCTF® 135 and NCTF® 135 HA between 2 and 25°C (36-77°F) in a dry place protected from light and moisture. Handle with caution.

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